


Protection-class BM-003 Respirator FFP2 is made for environments in which harmful and mutagenic particles can be found in the air we breathe. Respirator masks of this class provide a barrier to at least 95% of particles measuring larger than 0.3 micron and may be used in environments where hazardous substances exceed the OEL by up to a maximum of tenfold concentration. Our filter system, with a two layers of meltblown cloth, protects wearers from unpleasant odours and particles on top of the required breathing protection.



**Standard:** PPE marking, FFP2NR, CE, EN 149:2001+A1:2009, COVID-19 PPE-R/02.075-V2

**Filtration efficiency:** Above 95% (Equivalent to US N95 and EU FFP2 standard)

**Shelf life:** Two years

**Designation:** NR - Non reusable (single shift use only)

**Disposal:** Used products should be disposed as hazardous waste in accordance with national regulations.

**Materials:**

50gsm Spunbond Nonwoven

25gsm Meltblown Nonwoven

40-50gsm Hot air Nonwoven

25gsm Meltblown Nonwoven

25gsm Spunbond Nonwoven


**Mask weight:** 6 grams

# Technical Data sheet

Baltic Masks, BM – 003 respirator, FFP2 (v1.04)

## Packaging

2 pcs in a polybag, 10 pcs in a inner box, 48 inner boxes in a carton box



1 Size: 150mm (y) x 64mm (x) x 120mm (z)  
Weight: 110g  
Quantity: 1box ( 10pcs.)

2 Size: 490mm (y) x 400mm (x) x 320mm (z)  
Weight: 5.3kg  
Quantity: 1/48box ( 480pcs.)

3 Size: 980mm (y) x 800mm (x) x 1920mm (z)  
Weight: 128kg  
Quantity: 24/48box ( 11.520pcs.)

## Container loading

Container size: 20'DC

Loading: Carton box, free loading

Quantity: 499 carton / 239 520 pcs.

Container size: 40'DC


Loading: Carton box, free loading

Quantity: 1032 carton / 495 360 pcs.

Container size: 40'HC

Loading: Carton box, free loading

Quantity: 1110 carton / 532 800 pcs.



**Storage Conditions and Shelf Life** The filter masks have to be stored in their packets to retain their properties. Shelf life of unopened product is 2 years from date of manufacture when stored within temperature range of -20°C to +40°C and at less than 80% relative humidity.

12/5/2020

 **Usage instruction**

1. Before opening and touching the mask, wash your hands thoroughly with soap and water or alcohol/based hand sanitizer (>60% alcohol content).
2. Open the mask with both hands to the two sides.
3. Fit the mask on your face covering the mouth and nose. The straps should not be between the mask and your face.
4. Pull the straps above your ears. Adjust the mask until it fits comfortably.
5. Adjust the nose bridge stripe, until it completely fits to your nose bridge.
6. Cover the mask with both hands until the mask fits tightly around the face.

 **Precautions**

1. This product is not suitable for protection against harmful gases, steam and oily aerosols.
2. The product does not generate oxygen and is not suitable for use in anoxic environments, underwater operations, during evacuation or firefighting environments.
3. It is not recommended to use this protective mask if you do not feel a good adhesion between the face and the edge of the mask during use.
4. If found during use that the mask is damaged or the breathing resistance is too large, replace with a new one.
5. This product is for one-time use, masks should not be reused. Recommended use time is 4-8 hours.

*This filtering half mask is manufactured for COVID-protection only. This filtering half mask is not a PPE device for general use and shall not be used for purposes other than protection against COVID-19. This product shall not be used for more than one shift.*



Changxing Zhiwei Clothing Materials Co.,LTD Xingxin Industrial Park,  
Changxing county,Huzhou City,Zhejiang Province,China

12/5/2020